Bürkle GmbH

Rheinauenstr. 5 | 79415 Bad Bellingen
Telefon +49 7635 6121-0
info@buerkle.de
www.buerkle.de

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General safety requirements

- ▶ Read these directions for use carefully before using.
- ▶ Keep these directions for use where they will be accessible for all users.
- ▶ Failure to comply with the directions for use shall render all warranty conditions null and void and may lead to personal injury and material damage.
- ▶ The instrument must only be used by qualified skilled employees.

Safety instructions before using

- ▶ Follow the manufacturer's instructions and safety instructions of the medium you are sampling.
- ▶ When sampling hazardous substances, follow the technical rules on hazardous substances (in accordance with TRGS – (German) Technical Rules on Hazardous Substances) and wear appropriate protective equipment.
- ▶ The instrument is only suitable to take samples and not for continuous operation.
- ▶ Only use the instrument under adequate weather conditions.

WARNING!

Damage to the respiratory tract and other organs due to escaping or outgassing media.

- ▶ When sampling hazardous substances, follow the safety data sheet of the medium to be sampled.
- ▶ Wear the appropriate protective clothing.
- ▶ Do not use the instrument in frozen watercourses.
- ▶ Operate the instrument only with the human power of a person.
- ▶ Secure the sampling area before sampling (pay attention to humans, objects, ships, moving parts etc.).
- ▶ The instrument is not used for removing bulk goods.
- ▶ Do not use any force when operating the instrument.
- ▶ The instrument must only be used within the confines of applicable material resistance and for its designated purpose. Any other type of use will affect the warranty/guarantee.

NOTE!

- Before using, the instrument has to be cleaned to the extent of your (validated) cleaning process.
- ▶ Measures of manual cleaning have to be carried out in the cleaning/desinfection area and under consideration of occupational health and safety.
 - ▶ Pay attention to the material resistance.

CAUTION!

- Possible loss of the sample container during use.
- ▶ Always pay attention on the sampling area.

Safety instructions during and after use

- ▶ The instrument is suitable at a temperature range from -20 to 60 °C.
- ▶ Pay attention to have sufficient space during use: Radius of the telescopic rod with attached tool + 50 cm.
- ▶ Only perform lifting, lowering and slight lateral movements during use.

NOTE!

- The instrument is suitable for use with the telescopic rods specified in the accessories.
- ▶ Do not use any other tools than those specified in the accessories.

Layout and function

- ▶ The instrument is made of stainless steel V2A (1.4301).
- ▶ The instrument is suitable to take samples of liquids, low-viscosity media and pourable powders and slurries.
- ▶ The instrument is maintenance-free and suitable for steam sterilisation (up to 121°C) or sterilisation by flaming.

Parts in contact with medium

- ▶ Stainless steel V2A (1.4301)
- ▶ Aluminium (telescopic rod), for immersion depth > 50 cm
- ▶ Galvanised steel (snap-in joint), for immersion depth > 50 cm

Assembly of the telescopic rod

1. Connect the beaker to the telescopic rod using the snap-in joint ①.

CAUTION!

- Moving parts on the telescopic rod. Risk of crushing.
- ▶ Carefully operate the snap connection when connecting the tool.

2. Loosen screw connections of telescopic rod.
3. Extend the telescopic rod to the desired length.
4. Tighten the screw connections of telescopic rod by hand.

Safe sampling

WARNING!

Risk of slipping and falling due to wet or smooth floor.

- ▶ Ensure to have enough support and stability.

WARNING!

Breakage during sampling.

- ▶ Do not overload the beaker.
- ▶ Perform a visual check for material fatigue (cracks etc.) before sampling.
- ▶ Ensure that the tool is attached correctly before use (the snap-in joint has to be locked in both holes).

CAUTION!

Deformation or breakage of the telescopic rod.

- ▶ Avoid excessive lever movements.

CAUTION!

Bodily injury due to uncontrolled swinging around of the instrument and the telescopic rod.

- ▶ Pay attention to have sufficient space during use: Radius of the telescopic rod with attached tool + 50 cm.
- ▶ Disassemble the instrument for transport.

1. Take the sample using the beaker.
2. Empty the sample into a suitable collecting container.
3. Remove the beaker from the telescopic rod.
4. Clean the instrument.

Accessories

Type	Specification	Number	Item No.
Telescopic rod	65 - 120 cm	1	5355-0120
Telescopic rod	70 - 180 cm	1	5355-0180
Telescopic rod	125 - 250 cm	1	5355-0250
Telescopic rod	95 - 280 cm	1	5355-0280
Telescopic rod	115 - 300 cm	1	5355-0300
Telescopic rod	165 - 450 cm	1	5355-0450
Telescopic rod	175 - 600 cm	1	5355-0600
Scoop angular beaker	600 ml	1	5354-0050
Scoop angular beaker	1000 ml	1	5354-0100
Scoop angular beaker	2000 ml	1	5354-0200
Pendulum beaker PP	600 ml	1	5354-5050
Pendulum beaker PP	1000 ml	1	5354-5100
Pendulum beaker PP	2000 ml	1	5354-5200
Pendulum beaker AISI 304	1000 ml	1	5354-3000
Bottle holder	750 ml	1	5354-4100
Surface catch net	310 x 280 mm	1	5354-5020
Depth catch net	190 x 430 mm	1	5354-5010

More accessories can be found at www.buerkle.de

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