

Declaration of Conformity

on plastic materials and articles intended to come into contact with food

The on page 3 mentioned products made of LDPE are in compliance with the legal requirements of Regulation (EU) No. 10/2011, (incl. Regulation (EU) 2017/752), 1282/2011 and 1183/2012, Regulation (EU) 1935/2004 and Regulation (EC) 2023/2006 as amended. The total migration, as well as the specific migrations are below the legal limit value, according to the Regulation (EU) 10/2011.

In case the products contain restricted substances (SML / OM), it will comply with the limit values that are listed in Regulation (EU) 10/2011. For substances, that are not listed in the Union List (Annex 1 of Regulation (EU) 10/2011), apply the national provisions such as the recommendations of the BfR (the German federal government; Federal Institute for Risk Assessment).

The mentioned products can be used for long term storage of all foodstuff at room temperature or below and can be heated up to 70 °C for up to 2 hours or to 100 °C for up to 15 minutes.

Thus, on page 3 mentioned products with the above mentioned restrictions also comply with the provisions of Regulation (EU) No. 1935/2004 and the Food and Feed Code (LFGB).

Ratio of the area in contact with food to the volume by which the conformity of the material or article has been established:

6 dm²/kg

Functional barriers were not used.

The user has to convince himself of the suitability of the product for the intended filling material which goes beyond the requirements of the directives.

The polymer grade does contain one or more additive(s) that is/are subject to a restriction in food as referred to in Article 11.3 of EU Regulation 10/2011.

Note

For information purpose only

This note contains information relative to the presence of additives subject to a restriction according to Regulation (EU) 10/2011, as described in this statement.

Additive : Octadecyl 3-(3,5-di-tert-butyl-4-hydroxyphenyl)propionate

EC Ref. No: 68320

Max. conc.*: 760 ppm

Restriction: SML = 6 mg/kg food - Lipophilic substance

Additive: Polyethylene glycol

EC Ref. No: 76960

Max. conc.*: 300 ppm

Restriction: Dual use additive

* This information is provided for general guidance purposes only and the raw material supplier provides no guarantees or warranties in respect of this information and has no responsibility or liability for any use by any third party of this information.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of the certificate is 2284-08.

This document has been translated for informative purpose. In case of doubt the German version is valid.

19. November 2019



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



Article	Description
0300-3018	Dropper cap thread dia. 18 mm
0300-3025	Dropper cap thread dia. 25 mm
0300-3028	Dropper cap thread dia. 28 mm
0302-0010	Narrow-necked bottles LDPE 10 ml
0302-0020	Narrow-necked bottles LDPE 20 ml
0302-0030	Narrow-necked bottles LDPE 30 ml
0302-0050	Narrow-necked bottles LDPE 50 ml
0302-0100	Narrow-necked bottles LDPE 100 ml
0302-0250	Narrow-necked bottles LDPE 250 ml
0302-0500	Narrow-necked bottles LDPE 500 ml
0302-1000	Narrow-necked bottles LDPE 1000 ml
0302-2000	Narrow-necked bottles LDPE 2000 ml
0306-0010	Dropping bottles LDPE 10 ml
0306-0020	Dropping bottles LDPE 20 ml
0306-0030	Dropping bottles LDPE 30 ml
0306-0050	Dropping bottles LDPE 50 ml
0306-0100	Dropping bottles LDPE 100 ml
0306-0250	Dropping bottles LDPE 250 ml
0306-0500	Dropping bottles LDPE 500 ml
0306-1000	Dropping bottles LDPE 1000 ml
0310-0010	Wash bottle LDPE 100 ml
0310-0025	Wash bottle LDPE 250 ml
0310-0050	Wash bottle LDPE 500 ml
0310-0100	Wash bottle LDPE 1000 ml
0310-0200	Wash bottle LDPE 2000 ml
0314-0100	Tall-form bottles LDPE 100 ml
0314-0250	Tall-form bottles LDPE 250 ml
0314-0500	Tall-form bottles LDPE 500 ml
0314-1000	Tall-form bottles LDPE 1000 ml
0318-0050	Wide-necked bottle LDPE 50 ml
0318-0100	Wide-necked bottle LDPE 100 ml
0318-0250	Wide-necked bottle LDPE 250 ml
0318-0500	Wide-necked bottle LDPE 500 ml
0318-1000	Wide-necked bottle LDPE 1000 ml
0318-2000	Wide-necked bottle LDPE 2000 ml
0320-0250	Wide-necked bottle LDPE, brown 250 ml
0320-0500	Wide-necked bottle LDPE, brown 500 ml
0320-1000	Wide-necked bottle LDPE, brown 1000 ml